

## HUMAN CLDN7 FULL LENGTH PROTEIN

货号: 11061

产品全名: 人 CLDN7 Full Length 蛋白

规格: 10/50/100 µg

基因符号 CEPTL2; claudin-1; CLDN-7; CPETRL2; Hs.84359

目标蛋白: CLDN7

**UNIPROT ID:** O95471

**描述:** 人 CLDN7 full length protein-synthetic nanodisc

**背景:** This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. Differential expression of this gene has been observed in different types of malignancies, including breast cancer, ovarian cancer, hepatocellular carcinomas, urinary tumors, prostate cancer, lung cancer, head and neck cancers, thyroid carcinomas, etc.. Alternatively spliced transcript variants encoding different isoforms have been found.[provided by RefSeq, May 2010]

**物种/宿主:** HEK293

**分子量:** The human full length CLDN7 protein has a MW of 22.9 kDa

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.