

## CLEC2D (DMC394) IGG1 CHIMERIC MAB

货号: 28212

产品全名: CLEC2D(DMC394) IgG1 Chimeric 单抗

基因符号 C-type lectin domain family 2 member D;Lectin-like NK cell receptor;LLT-1;Osteoclast inhibitory lectin

描述: CLEC2D antibody(DMC394) IgG1 Chimeric 单抗

背景: This gene encodes a member of the natural killer cell receptor C-type lectin family. The encoded protein inhibits osteoclast formation and contains a transmembrane domain near the N-terminus as well as the C-type lectin-like extracellular domain. Several alternatively spliced transcript variants have been identified for this gene.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit:Human Fc chimeric IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 CLEC2D

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-CLEC2D (DMC394) on Expi293 cells transfected with human CLEC2D (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).