

HUMAN CX3CR1 FULL LENGTH PROTEIN

货号: 11096

产品全名: 人 CX3CR1 Full Length 蛋白

规格: 10/50/100 µg

基因符号 CCRL1; CMKBRL1; CMKDRI; GPR13; GPRV28; V28

目标蛋白: CX3CR1

UNIPROT ID: P49238

描述: 人 CX3CR1 full length protein-synthetic nanodisc

背景: Fractalkine is a transmembrane protein and chemokine involved in the adhesion and migration of leukocytes. The protein encoded by this gene is a receptor for fractalkine. The encoded protein also is a coreceptor for HIV-1, and some variations in this gene lead to increased susceptibility to HIV-1 infection and rapid progression to AIDS.

物种/宿主: HEK293

分子量: The human full length CX3CR1 protein has a MW of 40.4 kDa

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.