

CYSTATIN C MAB

产品名称: Cystatin C 小鼠单抗

货号: 26085

基因符号: Cystatin C

描述: Cystatin C 小鼠单抗

背景: Cystatin C is mainly used as a biomarker of kidney function. It is a potent inhibitor of cysteine proteases.

免疫原: Recombinant full-length Cystatin C protein, human origin.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:100-1:5000

WB: 1:100-1:1000

IHC: 1:100-1:1000

Concentration: 1 mg/ml

种属反应性: Mouse

形式: liquid

克隆性: Monoclonal

亚型: IgG

纯化: Purified from ascites


成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol
Species

种属反应性: Cystatin C antibody recognizes Cystatin C proteins of vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing.

免疫沉淀/蛋白质印迹:

Cystatin C protein



WB: Anti-Cystatin C mAb

Western blot analysis of recombinant full-length Cystatin C protein.

Purified His-tagged Cystatin C protein (货号:10157) was blotted with anti Cystatin C monoclonal antibody (货号: 26085) at 1:1000 dilution.