

**CYSTATIN C PAB**

产品名称: Cystatin C 兔多抗

货号: 21108

基因符号: Cystatin C

描述: Cystatin C 兔多抗

背景: Cystatin C, belongs to the cystatin family, is an inhibitor of cysteine proteinases, and is thought to serve an important physiological role as a local regulator of this enzyme activity. It is secreted. Cystatin C is expressed in the highest levels in the epididymis, vas deferens, brain, thymus, and ovary and the lowest in the submandibular gland. It is a potential cerebrospinal fluid marker for the diagnosis of Creutzfeldt-Jakob disease.

免疫原: A synthetic peptide from the internal region of Cystatin C, human origin.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:100-1:5000

WB: 1:100-1:1000

IHC: 1:100-1:1000

**Concentration:** 0.1 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG

纯化: Purified from serum


成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Cystatin C antibody recognizes Cystatin C proteins of vertebrates.

储存条件: Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing

蛋白质印迹:

**Cystatin C protein**



**WB: Anti-Cystatin C pAb**

**Western blot analysis of recombinant full length Cystatin C protein.**

Purified His-tagged Cystatin C protein(货号:10157) with anti-Cystatin C polyclonal anti.