

EMCN (DMC473) IGG1 CHIMERIC MAB

货号: 28249

产品全名: EMCN(DMC473) IgG1 Chimeric 单抗

基因符号 Endomucin; MUC-14

描述: EMCN antibody(DMC473) IgG1 Chimeric 单抗

背景: EMCN is a mucin-like sialoglycoprotein that interferes with the assembly of focal adhesion complexes and inhibits interaction between cells and the extracellular matrix (Kinoshita et al.; 2001 [PubMed 11418125]).[supplied by OMIM; Mar 2008]

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit:Human Fc chimeric IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 EMCN

成分: Lyophilized from sterile PBS, pH 7.4. 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-EMCN (DMC473) on Expi293 cells transfected with human EMCN (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
