

CGRPR MAB (ERENUMAB)

货号: 28157

产品全名: CGRPR(Erenumab) 单抗

基因符号 CALCRL

描述: CGRPR(Erenumab) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Homo sapiens

亚型: IgG2

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 CGRPR

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with $1\mu\text{g}/\text{mL}$ Anti-CGRPR(Erenumab biosimilar) mAb (28157) on Expi293 cells transfected with Human CGRPR&RAMPI protein (Blue histogram) and Expi293 transfected with Human CGRPR protein (Red histogram) or Human RAMPI protein(Orange histogram).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
