

EYA2 PAB

产品名称: EYA2

货号: 21128

基因符号: EYA2

描述: EYA2 兔多抗

背景: EYA2 (eyes absent homolog 2) functions as a tyrosine phosphatase that specifically dephosphorylates Tyr-142 of histone H2AX. EYA2 promotes efficient DNA repair and coactivates SIX1. Together with SIX1 and DACH2, EYA2 seems to be involved in myogenesis.

免疫原: A synthetic peptide from the internal region of human EYA2.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:1000-1:5000

WB: 1:100-1:1000

IHC: 1:100-1:1000

Concentration: 0.1 mg/ml

种属反应性: Rabbit

形式:

克隆性: Polyclonal

亚型: IgG

纯化: Purified from serum


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 59% glycerol
Species

种属反应性: EYA2 antibody recognizes EYA2 of vertebrates.

储存条件: Store at -20°C. Avoid repeated freezing and thawing

免疫沉淀/蛋白质印迹:

EYA2 protein



WB: Anti-EYA2 pAb

Western blot analysis of recombinant full length EYA2 protein.

Purified His-tagged EYA2 protein (货号:10195) was blotted with anti-EYA2 兔多抗 (Cat # 21128).