

## CD19 MAB (FMC63 MAB)

货号: 28095

产品全名: CD19(FMC63) 单抗

基因符号 CD19;4;VID3;GC12802

描述: CD19(FMC63) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Mus musculus

亚型: IgG2a

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Mouse CD19

成分: Lyophilized from sterile PBS, pH 7.4. 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1  $\mu\text{g}/\text{mL}$  Anti-CD19 mAb (28095) on Expi293 cells transfected with Human CD19 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---