

GA12 PAB

产品名称: $G\alpha_{12}$ 兔多抗

货号: 21013

基因符号: Gna12

描述: $G\alpha_{12}$ 兔多抗

背景: Heterotrimeric G proteins are essential cellular signal transducers. $G\alpha_{12}$ is a member of the $G_{12/13}$ subfamily of G proteins. Biochemical and physiological functions of $G\alpha_{12}$ are being explored.

免疫原: Recombinant full length $G\alpha_{12}$ 蛋白

经过测试的应用: ELISA, WB, IHC

推荐稀释比: ELISA: 1:1000-1:5000

WB: 1:100-1:2000

Concentration: 1 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG


纯化: Purified from serum

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: $G\alpha_{12}$ 兔多抗 recognizes $G\alpha_{12}$ of vertebrates.

储存条件: Store at -20°C . Avoid repeated freezing and thawing

蛋白质印迹:



WB: anti- $G\alpha_{12}$ pAb

Western blot analysis of $G\alpha_{12}$ proteins.

MEF whole-cell lysates were used for western blot with Anti- $G\alpha_{12}$ 兔多抗 (cat # 21013).

