

## GA13 PAB

产品名称:  $G\alpha_{13}$  兔多抗

货号: 21020

基因符号: GNA13

描述:  $G\alpha_{13}$  兔多抗

背景: Heterotrimeric G proteins are essential cellular signal transducers.  $G\alpha_{13}$  is one of the G proteins that could mediate cell migration and angiogenesis. Other biochemical and physiological functions of  $G\alpha_{13}$  are being explored.

免疫原: A synthetic peptide from a unique region of  $G\alpha_{13}$

经过测试的应用: ELISA, WB, IHC

推荐稀释比:

ELISA: 1:100-1:5000

WB: 1:100-1:2000

IHC: 1:50-1:100

**Concentration:** 0.2 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal


纯化: Purified from serum

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性:  $G\alpha_{13}$  兔多抗 recognizes  $G\alpha_{13}$  of vertebrates.

储存条件: Store at  $-20^{\circ}\text{C}$ . Avoid repeated freezing and thawing

蛋白质印迹:



**WB: anti-Gα13 pAb**

**Western blot analysis of Gα13 proteins.**

Purified Gα13 protein was blot with Anti-Gα13 兔多抗 (Cat # 21020).