

GA14 MAB

产品名称: $G\alpha_{14}$

货号: 26027

基因符号: GNAI4

描述: $G\alpha_{14}$ 小鼠单抗

背景: Heterotrimeric G proteins are essential cellular signal transducers. $G\alpha_{14}$ is a member of the Gq subfamily of G proteins. $G\alpha_{14}$ is known to stimulate phospholipase C. Other biochemical and physiological functions of $G\alpha_{14}$ are being explored.

免疫原: A synthetic peptide from the C-terminal region of human $G\alpha_{14}$

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比: ELISA 1:100–1:5000

WB 1:100–1:2000

IHC 1:50–1:100

Concentration: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG


纯化: Purified from ascites

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: $G\alpha_{14}$ antibody recognizes $G\alpha_{14}$ of vertebrates.

储存条件: Store at -20°C . Avoid repeated freezing and thawing

蛋白质印迹:



WB: anti-Gna14 mAb

HepG2 whole-cell lysate was used for western blot (cat# 26027).