

HTRA1 MAB (GALEGENIMAB)

货号: 28141

产品全名: HTRA1(Galegenimab) 单抗

基因符号 L56;Serine protease II

描述: HTRA1(galegenimab) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 HTRA1

成分: Lyophilized from sterile PBS, pH 7.4. 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 15 $\mu\text{g}/\text{mL}$ Anti-HTRA1(galegenimab biosimilar) mAb (28141) on Expi293 cells transfected with Human HTRA1 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
