

## GA<sub>o</sub> PAB

产品名称: GA<sub>o</sub> 兔多抗

货号: 21015

基因符号: Gnao

描述: GA<sub>o</sub> 兔多抗

背景: Heterotrimeric G proteins are critical cellular signal transducers. GA<sub>o</sub> belongs to the Gi sub-family of G proteins that could mediate the inhibition of adenylyl cyclases. Other biochemical and physiological functions of GA<sub>o</sub> proteins are being explored.

免疫原: Recombinant full length protein of GA<sub>o</sub>.

经过测试的应用: ELISA, WB, IHC

推荐稀释比: ELISA 1:1000-1:5000

WB 1:100-1:2000

**Concentration:** 1 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG


纯化: Purified from serum

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: GA<sub>o</sub> antibody recognizes GA<sub>o</sub> of vertebrates.

储存条件: Store at -20°C. Avoid repeated freezing and thawing

蛋白质印迹:



WB: anti-Gα<sub>o</sub> pAb

### Western blot analysis of Go proteins.

Whole-cell lysate of MEF cells or purified Gα<sub>o</sub> protein were used for western blots with the polyclonal anti Gα<sub>o</sub> antibody (Cat # 21015).