

## GAQ PAB

产品名称: Ga<sub>q</sub> 兔多抗

货号: 21014

基因符号: Gnaq

描述: Gaq 兔多抗

背景: Heterotrimeric G proteins are essential cellular signal transducers. Gaq is a member of the Gq subfamily of G proteins. Gaq is known to stimulate phospholipase C. Other biochemical and physiological functions of Gaq are being explored.

免疫原: Recombinant full length Gaq 蛋白

经过测试的应用: ELISA, WB, IHC

推荐稀释比: ELISA 1:1000-1:5000 WB 1:100-1:2000

**Concentration:** 1 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal


亚型: IgG

纯化: Purified from serum

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Gaq 兔多抗 recognizes Gaq of vertebrates.

储存条件: Store at -20°C. Avoid repeated freezing and thawing



蛋白质印迹: WB: anti-Ga<sub>q</sub> pAb **Western blot analysis of Gaq proteins.** MEF whole-cell lysates were used for western blot with Anti-Gaq 兔多抗 (cat # 21014).