

GAS MAB

产品名称: $G\alpha_s$ 小鼠单抗

货号: 26006

基因符号: Gnas

描述: $G\alpha_s$ 小鼠单抗

背景: Heterotrimeric G proteins are critical cellular signal transducers. $G\alpha_s$ represents one sub-family of G proteins that could mediate the stimulation of adenylyl cyclases.

免疫原: Recombinant full length protein of $G\alpha_s$

经过测试的应用: ELISA, WB, IHC

推荐稀释比: ELISA 1:1000-1:5000, WB 1:1000-1:2000, IHC 1:100-1:200

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG


纯化: Purified from ascites

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 50% glycerol

种属反应性: $G\alpha_s$ antibody recognizes $G\alpha_s$ of vertebrates

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

蛋白质印迹:



Western blot analysis of $G\alpha_s$ proteins. MEF whole-cell lysate was used for western blots.