

**GAS PAB**

产品名称:  $G\alpha_s$  兔多抗

货号: 21007

基因符号: Gnas

描述:  $G\alpha_s$  兔多抗

背景: Heterotrimeric G proteins are critical cellular signal transducers.  $G\alpha_s$  represents one sub-family of G proteins that could mediate the stimulation of adenylyl cyclases.

免疫原: Recombinant full length protein of  $G\alpha_s$

经过测试的应用: ELISA, WB, IHC

推荐稀释比: ELISA 1:100-1:5000 WB 1:100-1:2000

**Concentration:** 1 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG


纯化: Purified from serum

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性:  $G\alpha_s$  antibody recognizes  $G\alpha_s$  of vertebrates

储存条件: Store at  $-20^{\circ}\text{C}$ . Avoid repeated freezing and thawing

蛋白质印迹:



WB: anti- $G\alpha_s$  pAb

**Western blot analysis of  $G\alpha_s$  proteins.**

MEF whole-cell lysate was used for western blots.

