

GAZ MAB

产品名称: $G\alpha_z$ 小鼠单抗

货号: 26011

基因符号: GNAZ

描述: $G\alpha_z$ 小鼠单抗

背景: Heterotrimeric G proteins are critical cellular signal transducers. $G\alpha_z$ belongs to the G_i sub-family of G proteins that could mediate the inhibition of adenylyl cyclases. Other biochemical and physiological functions of $G\alpha_z$ proteins are being explored.

免疫原: Recombinant full length protein of $G\alpha_z$.

经过测试的应用: ELISA, WB, IHC

推荐稀释比: ELISA 1:100-1:5000; WB 1:100-1:2000

Concentration: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG


纯化: Purified from ascites

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: $G\alpha_z$ antibody recognizes $G\alpha_z$ of vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

蛋白质印迹:



WB: anti- $G\alpha_z$ mAb

Western blot analysis of $G\alpha_z$ proteins. HEK-293 whole-cell lysate was used for western blots.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
