

## HUMAN GLP2R FULL LENGTH PROTEIN

货号: 11094

产品全名: 人 GLP2R Full Length 蛋白

规格: 10/50/100 µg

基因符号 GLP-2-R; GLP-2R

目标蛋白: GLP2R

**UNIPROT ID:** O95838

**描述:** 人 GLP2R full length protein-synthetic nanodisc

**背景:** A G protein-coupled receptor that is closely related to the glucagon receptor and binds to glucagon-like peptide-2 (GLP2). Signalling through GLP2 stimulates intestinal growth and increases villus height in the small intestine, concomitant with increased crypt cell proliferation and decreased enterocyte apoptosis.

**物种/宿主:** HEK293

**分子量:** The human full length GLP2R protein has a MW of 63.0 kDa

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.