

HUMAN GPR65 FULL LENGTH PROTEIN

货号: 11084

产品全名: 人 GPR65 Full Length 蛋白

规格: 10/50/100 µg

基因符号 hTDAG8; TDAG8

目标蛋白: GPR65

UNIPROT ID: Q8IYL9

描述: 人 GPR65 full length protein-synthetic nanodisc

背景: Receptor for the glycosphingolipid psychosine (PSY) and several related glycosphingolipids (PubMed:11309421). Plays a role in immune response by maintaining lysosome function and supporting phagocytosis-mediated intracellular bacteria clearance (PubMed:27287411). May have a role in activation-induced cell death or differentiation of T-cells (By similarity).[UniProtKB/Swiss-Prot Function]

物种/宿主: HEK293

分子量: The human full length GPR65 protein has a MW of 39.3 kDa

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.