

GPR87 (DMC478) IGG1 CHIMERIC MAB

货号: 28254

产品全名: GPR87(DMC478) IgG1 Chimeric 单抗

基因符号 FKSG78; GPR95; KPG_002

描述: GPR87 antibody(DMC478) IgG1 Chimeric 单抗

背景: This gene encodes a G protein-coupled receptor and is located in a cluster of G protein-couple receptor genes on chromosome 3. The encoded protein has been shown to be overexpressed in lung squamous cell carcinoma (PMID:18057535) and regulated by p53 (PMID:19602589). [provided by RefSeq; Nov 2011]

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit:Human Fc chimeric IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 GPR87

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-GPR87 (DMC478) on Expi293 cells transfected with human GPR87 (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram)



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
