

HUMAN GPRC5D FULL LENGTH PROTEIN

货号: 11002

产品全名: 人 GPRC5D Full Length 蛋白

规格: 10/50/100 µg

基因符号 GPRC5D

目标蛋白: GPRC5D

UNIPROT ID: Q9NZD1

描述: 人 GPRC5D full length protein-synthetic nanodisc

背景: The protein encoded by this gene is a member of the G protein-coupled receptor family. Recent studies demonstrate that GPRC5D is expressed on malignant bone marrow plasma cells, whereas normal tissue expression is limited to the hair follicle. It may represent a potential target for effector-cell-mediated therapy to treat plasma-cell disorders like MM.

物种/宿主: HEK293

分子量: The human full length GPRC5D Protein has a MW of 38.6 kDa

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure1. Elisa plates were added with Flag Tag GPRC5D-Nanodisc (0.5 µg/per well) on an anti-Flag monoclonal antibody pre-coated (0.5 µg/per well) plate. Serial diluted anti-GPRC5D monoclonal antibody (28421) solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for anti-GPRC5D monoclonal antibody binding with GPRC5D-Nanodisc is 32.86ng/ml.




Figure2. Human GPRC5D-Nanodisc, Flag Tag on SDS-PAGE