

GUCY2C (DMC387) IGG1 CHIMERIC MAB

货号: 28205

产品全名: GUCY2C(DMC387) IgG1 Chimeric 单抗

基因符号 DIAR6; GC-C; GUC2C; MECIL; MUCIL; STAR

描述: GUCY2C antibody(DMC387) IgG1 Chimeric 单抗

背景: This gene encodes a transmembrane protein that functions as a receptor for endogenous peptides guanylin and uroguanylin; and the heat-stable E. coli enterotoxin. The encoded protein activates the cystic fibrosis transmembrane conductance regulator. Mutations in this gene are associated with familial diarrhea (autosomal dominant) and meconium ileus (autosomal recessive).

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit:Human Fc chimeric IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 GUCY2C

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-GUCY2C (DMC387) on Expi293 cells transfected with human GUCY2C (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).