

HUMAN APOH PROTEIN

货号: 12270

产品全名: 人 APOH 蛋白

规格: 10/50/100 µg

基因符号 BG;B2G1;B2GPI

目标蛋白: APOH

UNIPROT ID: P02749

描述: Recombinant human APOH Protein with C-terminal 6xHis tag

背景: Apolipoprotein H, also known as beta-2-glycoprotein I, is a component of circulating plasma lipoproteins. It has been implicated in a variety of physiologic pathways including lipoprotein metabolism, coagulation, hemostasis, and the production of antiphospholipid autoantibodies. APOH may be a required cofactor for anionic phospholipid binding by the antiphospholipid autoantibodies found in sera of many patients with lupus and primary antiphospholipid syndrome (APS). The anti-beta (2) glycoprotein I antibodies from APS patients, mediate inhibition of activated protein C which has anticoagulant properties. Because beta-2-GPI is the main autoantigen in patients with APS, the disruption of this pathway by autoantibodies may be an important mechanism for thrombosis in patients with APS.[provided by RefSeq, Dec 2019]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 37.1 kDa after removal of the signal peptide. The apparent molecular mass of APOH-His is approximately 35-70 kDa due to glycosylation.

分子特征: APOH(Gly20-Cys345) 6×His tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human APOH Protein, His Tag on SDS-PAGE under reducing condition.