

HUMAN AXL PROTEIN, MFC TAG

货号: 11333

产品全名: 人 AXL 蛋白

规格: 10/50/100 µg

基因符号 AXL;FO

目标蛋白: AXL

UNIPROT ID: P30530

描述: Recombinant Human AXL with C-terminal mouse Fc tag

背景: The protein encoded by this gene is a member of the Tyro3-Axl-Mer (TAM) receptor tyrosine kinase subfamily. The encoded protein possesses an extracellular domain which is composed of two immunoglobulin-like motifs at the N-terminal, followed by two fibronectin type-III motifs. It transduces signals from the extracellular matrix into the cytoplasm by binding to the vitamin K-dependent protein growth arrest-specific 6 (Gas6). This gene may be involved in several cellular functions including growth, migration, aggregation and anti-inflammation in multiple cell types. Alternative splicing results in multiple transcript variants of this gene.

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 72.0 kDa after removal of the signal peptide. The apparent molecular mass of AXL-mFc is approximately 100-130 kDa due to glycosylation.

分子特征: AXL(Ala26-Pro449) mFc(Pro99-Lys330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human AXL Protein, mFc Tag on SDS-PAGE under reducing condition.