

HUMAN B7-H3 PROTEIN, MFC TAG

货号: 11543

产品全名: 人 B7-H3 蛋白

规格: 10/50/100 µg

基因符号 B7H3;CD276;4Ig-B7-H3

目标蛋白: B7-H3

UNIPROT ID: Q5ZPR3

描述: Recombinant human B7-H3 protein with C-terminal mouse Fc tag

背景: The protein encoded by this gene belongs to the immunoglobulin superfamily, and thought to participate in the regulation of T-cell-mediated immune response. Studies show that while the transcript of this gene is ubiquitously expressed in normal tissues and solid tumors, the protein is preferentially expressed only in tumor tissues. Additionally, it was observed that the 3' UTR of this transcript contains a target site for miR29 microRNA, and there is an inverse correlation between the expression of this protein and miR29 levels, suggesting regulation of expression of this gene product by miR29. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 72.8 kDa after removal of the signal peptide. The apparent molecular mass of B7-H3-mFc is approximately 70-100 kDa due to glycosylation.

分子特征: B7-H3(Gly27-Thr461) mFc(Pro99-Lys330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human B7-H3 Protein, mFc Tag on SDS-PAGE under reducing condition.