

HUMAN BCHE PROTEIN, HIS TAG

货号: 11872

产品全名: 人 BCHE 蛋白

规格: 10/50/100 µg

基因符号 BCHE;CHE1;CHE2;E1

目标蛋白: BCHE

UNIPROT ID: P06276

描述: Recombinant Human BCHE Protein with C-terminal 6xHis tag

背景: This gene encodes a cholinesterase enzyme and member of the type-B carboxylesterase/lipase family of proteins. The encoded enzyme exhibits broad substrate specificity and is involved in the detoxification of poisons including organophosphate nerve agents and pesticides, and the metabolism of drugs including cocaine, heroin and aspirin. Humans homozygous for certain mutations in this gene exhibit prolonged apnea after administration of the muscle relaxant succinylcholine. [provided by RefSeq, Jul 2016]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 65.9 kDa after removal of the signal peptide. The apparent molecular mass of BCHE-His is approximately 70-130 kDa due to glycosylation.

分子特征: BCHE(Glu29-Leu602) 6xHis tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human BCHE Protein, His Tag on SDS-PAGE under reducing condition.