

HUMAN CD2 PROTEIN, HFC TAG

货号: 11392

产品全名: 人 CD2 蛋白

规格: 10/50/100 µg

基因符号 LFA-2;SRBC;T11

目标蛋白: CD2

UNIPROT ID: P06729

描述: Recombinant Human CD2 Protein with C-terminal human Fc tag

背景: The protein encoded by this gene is a surface antigen found on all peripheral blood T-cells. The encoded protein interacts with LFA3 (CD58) on antigen presenting cells to optimize immune recognition. A locus control region (LCR) has been found in the 3' flanking sequence of this gene. [provided by RefSeq, Jun 2016]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 47.4 kDa after removal of the signal peptide. The apparent molecular mass of CD2-hFc is approximately 55-70 kDa due to glycosylation.

分子特征: CD2(Lys25-Asp209) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 90% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human CD2 Protein, His Tag on SDS-PAGE under reducing condition.