

HUMAN CD200R1 PROTEIN, HIS TAG

货号: 11608

产品全名: 人 CD200R1 蛋白

规格: 10/50/100 µg

基因符号 CD200R;HCRTR2;MOX2R;OX2R

目标蛋白: CD200R1

UNIPROT ID: Q8TD46

描述: Recombinant Human CD200R1 Protein with C-terminal 6xHis tag

背景: This gene encodes a receptor for the OX-2 membrane glycoprotein. Both the receptor and substrate are cell surface glycoproteins containing two immunoglobulin-like domains. This receptor is restricted to the surfaces of myeloid lineage cells and the receptor-substrate interaction may function as a myeloid downregulatory signal. Mouse studies of a related gene suggest that this interaction may control myeloid function in a tissue-specific manner. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jul 2008]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 25.0 kDa after removal of the signal peptide. The apparent molecular mass of CD200R1-His is approximately 35-70 kDa due to glycosylation.

分子特征: CD200R1(Met29-Leu243) 6×His tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human CD200R1 Protein, His Tag on SDS-PAGE under reducing condition.