

## HUMAN CD23 PROTEIN, HIS TAG

货号: 11440

产品全名: 人 CD23 蛋白

规格: 10/50/100 µg

基因符号 BLAST-2;CD23;CD23A;CLEC4J;FCE2;IGEBF

目标蛋白: CD23

**UNIPROT ID:** P06734

**描述:** Recombinant human CD23 protein with N-terminal 6xHis tag

**背景:** The protein encoded by this gene is a B-cell specific antigen, and a low-affinity receptor for IgE. It has essential roles in B cell growth and differentiation, and the regulation of IgE production. This protein also exists as a soluble secreted form, then functioning as a potent mitogenic growth factor. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.[provided by RefSeq, Jul 2011]

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 31.8 kDa after removal of the signal peptide. The apparent molecular mass of His-CD23 is approximately 35-55 kDa due to glycosylation.

**分子特征:** 6×His tag CD23(Asp48-Ser321)

**纯化:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human CD23 Protein, His Tag on SDS-PAGE under reducing condition.