

HUMAN CD235A PROTEIN, HFC TAG

货号: 11506

产品全名: 人 CD235A 蛋白

规格: 10/50/100 µg

基因符号 Glycophorin-A;CD235a;PAS-2;GYPA;GPA

目标蛋白: CD235A

UNIPROT ID: P02724

描述: Recombinant human CD235A protein with C-terminal human Fc tag

背景: Glycophorins A (GYPA) and B (GYPB) are major sialoglycoproteins of the human erythrocyte membrane which bear the antigenic determinants for the MN and Ss blood groups. In addition to the M or N and S or s antigens that commonly occur in all populations, about 40 related variant phenotypes have been identified. These variants include all the variants of the Miltenberger complex and several isoforms of Sta, as well as Dantu, Sat, He, Mg, and deletion variants Ena, S-s-U- and Mk. Most of the variants are the result of gene recombinations between GYPA and GYPB. [provided by RefSeq, Jul 2008]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 33.3 kDa after removal of the signal peptide. The apparent molecular mass of CD235A-hFc is approximately 35-55 kDa due to glycosylation.

分子特征: CD235A(His28-Glu91) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human CD235A Protein, hFc Tag on SDS-PAGE under reducing condition.