

## HUMAN CD26 PROTEIN, HFC TAG

货号: 11662

产品全名: 人 CD26 蛋白

规格: 10/50/100 µg

基因符号 ADABP;ADCP2;CD26;DPPIV;TP103

目标蛋白: ADABP

**UNIPROT ID:** P27487

**描述:** Recombinant human CD26 protein with C-terminal human Fc tag

**背景:** The DPP4 gene encodes dipeptidyl peptidase 4, which is identical to adenosine deaminase complexing protein-2, and to the T-cell activation antigen CD26. It is an intrinsic type II transmembrane glycoprotein and a serine exopeptidase that cleaves X-proline dipeptides from the N-terminus of polypeptides. Dipeptidyl peptidase 4 is highly involved in glucose and insulin metabolism, as well as in immune regulation. This protein was shown to be a functional receptor for Middle East respiratory syndrome coronavirus (MERS-CoV), and protein modeling suggests that it may play a similar role with SARS-CoV-2, the virus responsible for COVID-19. [provided by RefSeq, Apr 2020]

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 116.71 kDa after removal of the signal peptide.

**分子特征:** ADABP(Asn29-Pro766) hFc(Glu99-Ala330)

**纯化:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.