

HUMAN CD98 PROTEIN, HIS TAG

货号: 11963

产品全名: 人 CD98 蛋白

规格: 10/50/100 µg

基因符号 SLC3A2;F2;MDU1;4F2HC;4T2HC;NACAE;CD98HC

目标蛋白: CD98

UNIPROT ID: P08195

描述: Recombinant Human CD98 Protein with N-terminal 6XHis tag

背景: This gene is a member of the solute carrier family and encodes a cell surface, transmembrane protein. The protein exists as the heavy chain of a heterodimer, covalently bound through di-sulfide bonds to one of several possible light chains. The encoded transporter plays a role in regulation of intracellular calcium levels and transports L-type amino acids. Alternatively spliced transcript variants, encoding different isoforms, have been characterized. [provided by RefSeq, Nov 2010]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 47.7 kDa after removal of the signal peptide. The apparent molecular mass of His-CD98 is approximately 55-100 kDa due to glycosylation.

分子特征: 6×His tag CD98(Arg206-Ala630)

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human CD98 Protein, His Tag on SDS-PAGE under reducing condition.