

## HUMAN CGRP PROTEIN, HFC TAG

货号: 11962

产品全名: 人 CGRP 蛋白

规格: 10/50/100 µg

基因符号 CALCI;CALCA;CGRP-alpha;CGRP-I;CGRPI;CT;KC;PCT

目标蛋白: CGRP

**UNIPROT ID:** P06881

**描述:** Recombinant Human CGRP Protein with C-terminal human Fc tag

**背景:** This gene encodes the peptide hormones calcitonin, calcitonin gene-related peptide and katalcalcin by tissue-specific alternative RNA splicing of the gene transcripts and cleavage of inactive precursor proteins. Calcitonin is involved in calcium regulation and acts to regulate phosphorus metabolism. Calcitonin gene-related peptide functions as a vasodilator and as an antimicrobial peptide while katalcalcin is a calcium-lowering peptide. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Aug 2014]

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 29.9 kDa after removal of the signal peptide. The apparent molecular mass of CGRP-hFc is approximately 25-35 kDa due to glycosylation.

**分子特征:** CGRP(Ala83-Phe119) hFc(Glu99-Ala330)

**纯化:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human CGRP Protein, hFc Tag on SDS-PAGE under reducing condition.