

## HUMAN CXCR5 PROTEIN, HFC TAG

货号: 11616

产品全名: 人 CXCR5 蛋白

规格: 10/50/100 µg

基因符号: B1r1;CXC-R5;CXCR-5;Gpcr6;MDR15

目标蛋白: CXCR5

**UNIPROT ID:** P32302

**描述:** Recombinant Human CXCR5 with C-terminal human Fc tag

**背景:** Cytokine receptor that binds to B-lymphocyte chemoattractant (BLC). Involved in B-cell migration into B-cell follicles of spleen and Peyer patches but not into those of mesenteric or peripheral lymph nodes.[UniProtKB/Swiss-Prot Function]

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 32.3 kDa after removal of the signal peptide. The apparent molecular mass of CXCR5-hFc is approximately 35-55 kDa due to glycosylation.

**分子特征:** CXCR5(Met1-Phe53) hFc(Glu99-Ala330)

**纯化:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human CXCR4 Protein, hFc Tag on SDS-PAGE under reducing condition.