

HUMAN EPHA5 PROTEIN, HIS TAG

货号: 11411

产品全名: 人 EPHA5 蛋白

规格: 10/50/100 µg

基因符号 CEK7;EHK-1;EHK1;EK7;HEK7;TYRO4

目标蛋白: EPHA5

UNIPROT ID: P54756

描述: Recombinant human EPHA5 protein with C-terminal 6xHis tag

背景: This gene belongs to the ephrin receptor subfamily of the protein-tyrosine kinase family. EPH and EPH-related receptors have been implicated in mediating developmental events, particularly in the nervous system. Receptors in the EPH subfamily typically have a single kinase domain and an extracellular region containing a Cys-rich domain and 2 fibronectin type III repeats. The ephrin receptors are divided into 2 groups based on the similarity of their extracellular domain sequences and their affinities for binding ephrin-A and ephrin-B ligands. Alternatively spliced transcript variants encoding different isoforms have been described.

[provided by RefSeq, Aug 2013]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 61.6 kDa after removal of the signal peptide. The apparent molecular mass of EPHA5-His is approximately 70-100 kDa due to glycosylation.

分子特征: EPHA5(Pro25-Pro573) 6×His tag

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human EPHA5 Protein, His Tag on SDS-PAGE under reducing condition.