

HUMAN GPRC5C PROTEIN, MFC TAG

货号: 11776

产品全名: 人 GPRC5C 蛋白

规格: 10/50/100 µg

基因符号 RAIG-3

目标蛋白: GPRC5C

UNIPROT ID: Q9NQ84

描述: Recombinant human GPRC5C protein with C-terminal mouse Fc tag

背景: The protein encoded by this gene is a member of the type 3 G protein-coupled receptor family. Members of this superfamily are characterized by a signature 7-transmembrane domain motif. The specific function of this protein is unknown; however, this protein may mediate the cellular effects of retinoic acid on the G protein signal transduction cascade. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 28.8 kDa after removal of the signal peptide. The apparent molecular mass of GPRC5C-mFc is approximately 25-35 kDa due to glycosylation.

分子特征: GPRC5C(Val27-Gly50) mFc(Pro99-Lys330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human GPRC5C Protein, mFc Tag on SDS-PAGE under reducing condition.