

HUMAN IBSP PROTEIN, HIS TAG

货号: 11954

产品全名: 人 IBSP 蛋白

规格: 10/50/100 µg

基因符号 BSP;BNSP;SP-II;BSP-II

目标蛋白: IBSP

UNIPROT ID: P21815

描述: Recombinant Human IBSP Protein with C-terminal 6xHis tag

背景: The protein encoded by this gene is a major structural protein of the bone matrix. It constitutes approximately 12% of the noncollagenous proteins in human bone and is synthesized by skeletal-associated cell types, including hypertrophic chondrocytes, osteoblasts, osteocytes, and osteoclasts. The only extraskelatal site of its synthesis is the trophoblast. This protein binds to calcium and hydroxyapatite via its acidic amino acid clusters, and mediates cell attachment through an RGD sequence that recognizes the vitronectin receptor. [provided by RefSeq, Jul 2008]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 34.3 kDa after removal of the signal peptide. The apparent molecular mass of IBSP-His is approximately 55–70 kDa due to glycosylation.

分子特征: IBSP(Phe17–Gln317) 6×His tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at –20°C to –80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at –80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human IBSP Protein, His Tag on SDS-PAGE under reducing condition.