

## HUMAN IFN ALPHA6 (C-6HIS) PROTEIN

货号: 12040

产品全名: 人 IFN Alpha6 (C-6His) 蛋白

规格: 10/50/100 µg

基因符号 Interferon Alpha-6;IFN-Alpha-6;Interferon Alpha-54;Interferon Alpha-K;IeIF K;IFNA6

目标蛋白: IFN

**UNIPROT ID:** P05013

**描述:** Recombinant Human Interferon Alpha-6 is produced by our Mammalian expression system and the target gene encoding Ser21-Glu189 is expressed with a 6His tag at the C-terminus.

**背景:** Interferon α-6 (IFN-α6) is a secreted protein which belongs to the α/β interferon family. IFN-α6 is produced by macrophages, expressed at low level, only 1.0% of the average gene in this release. IFN-α6 contains interferon alpha, beta and delta domain. IFN-α has antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

物种/宿主: HEK293

分子量: 21.1 KDa

分子特征: Not available

纯化: Greater than 95% as determined by reducing SDS-PAGE.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Greater than 95% as determined by reducing SDS-PAGE.