

HUMAN IL23R PROTEIN, HIS TAG

货号: 11911

产品全名: 人 IL23R 蛋白

规格: 10/50/100 µg

基因符号 IL-23 receptor;IL-23R

目标蛋白: IL23R

UNIPROT ID: Q5VWK5

描述: Recombinant Human IL23R Protein with C-terminal 6xHis tag

背景: The protein encoded by this gene is a subunit of the receptor for IL23A/IL23. This protein pairs with the receptor molecule IL12RB1/IL12Rbeta1, and both are required for IL23A signaling. This protein associates constitutively with Janus kinase 2 (JAK2), and also binds to transcription activator STAT3 in a ligand-dependent manner. [provided by RefSeq, Jul 2008]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 37.4 kDa after removal of the signal peptide.

分子特征: IL23R(Val36-Asp353) 6×His tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human IL23R Protein, His Tag on SDS-PAGE under reducing condition.