

## HUMAN LRIG1 PROTEIN, HIS TAG

货号: 11523

产品全名: 人 LRIG1 蛋白

规格: 10/50/100 µg

基因符号 Leucine-rich repeats and immunoglobulin-like domains protein 1;IG-1

目标蛋白: LRIG1

**UNIPROT ID:** Q96JA1

**描述:** Recombinant human LRIG1 protein with C-terminal 6xHis tag

**背景:** Acts as a feedback negative regulator of signaling by receptor tyrosine kinases, through a mechanism that involves enhancement of receptor ubiquitination and accelerated intracellular degradation.[UniProtKB/Swiss-Prot Function]

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 84 kDa after removal of the signal peptide. The apparent molecular mass of LRIG1-His is approximately 100-130 kDa due to glycosylation.

**分子特征:** LRIG1(Ala35-Thr792) 6xHis tag

**纯化:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human LRIG1 Protein, His Tag on SDS-PAGE under reducing condition.