

HUMAN LY6H PROTEIN, HFC TAG

货号: 11908

产品全名: 人 LY6H 蛋白

规格: 10/50/100 µg

基因符号 NMLY6

目标蛋白: LY6H

UNIPROT ID: O94772

描述: Recombinant Human LY6H Protein with C-terminal human Fc tag

背景: Believed to act as a modulator of nicotinic acetylcholine receptors (nAChRs) activity. In vitro inhibits alpha-3:beta-4-containing nAChRs maximum response. May play a role in the intracellular trafficking of alpha-7-containing nAChRs and may inhibit their expression at the cell surface. Seems to inhibit alpha-7/CHRNA7 signaling in hippocampal neurons.[UniProtKB/Swiss-Prot Function]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 36.0 kDa after removal of the signal peptide. The apparent molecular mass of LY6H-hFc is approximately 35-55 kDa due to glycosylation.

分子特征: LY6H(Leu26-Gly115) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human LY6H Protein, hFc Tag on SDS-PAGE under reducing condition.