

HUMAN LYPD3(31-325) PROTEIN, HIS TAG

货号: 11891

产品全名: 人 LYPD3(31-325) 蛋白

规格: 10/50/100 µg

基因符号 MIG-C4;4.4A

目标蛋白: LYPD3

UNIPROT ID: O95274

描述: Recombinant Human LYPD3(31-325) Protein with C-terminal 6xHis tag

背景: Supports cell migration. May be involved in urothelial cell-matrix interactions. May be involved in tumor progression.[UniProtKB/Swiss-Prot Function]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 31.7 kDa after removal of the signal peptide. The apparent molecular mass of LYPD3(31-325)-His is approximately 35-130 kDa due to glycosylation.

分子特征: LYPD3(Leu31-Gly325) 6×His tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human LYPD3(31-325) Protein, His Tag on SDS-PAGE under reducing condition.