

HUMAN MFAP5 PROTEIN, HFC TAG

货号: 11582

产品全名: 人 MFAP5 蛋白

规格: 10/50/100 µg

基因符号 MFAP-5;MP25;MAGP-2;MAGP2

目标蛋白: MFAP5

UNIPROT ID: Q13361

描述: Recombinant human MFAP5 protein with C-terminal human Fc tag

背景: This gene encodes a 25-kD microfibril-associated glycoprotein which is a component of microfibrils of the extracellular matrix. The encoded protein promotes attachment of cells to microfibrils via alpha-V-beta-3 integrin. Deficiency of this gene in mice results in neutropenia. Alternate splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2014]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 43.4 kDa after removal of the signal peptide. The apparent molecular mass of MFAP5-hFc is approximately 35-55 kDa due to glycosylation.

分子特征: MFAP5(Ile22-Leu173) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human MFAP5 Protein, hFc Tag on SDS-PAGE under reducing condition.