

## HUMAN MMP9(20-469) PROTEIN, HIS TAG

货号: 11501

产品全名: 人 MMP9(20-469) 蛋白

规格: 10/50/100 µg

基因符号 GELB;CLG4B

目标蛋白: MMP9

**UNIPROT ID:** P14780

**描述:** Recombinant human MMP9 Protein with C-terminal 6xHis tag

**背景:** Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. The enzyme encoded by this gene degrades type IV and V collagens. Studies in rhesus monkeys suggest that the enzyme is involved in IL-8-induced mobilization of hematopoietic progenitor cells from bone marrow, and murine studies suggest a role in tumor-associated tissue remodeling. [provided by RefSeq, Jul 2008]

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 48.5 kDa after removal of the signal peptide. The apparent molecular mass of MMP9-His is approximately 55-70 kDa due to glycosylation.

**分子特征:** MMP9(Ala20-Pro469) 6×His tag

**纯化:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human MMP9, His Tag on SDS-PAGE under reducing condition.