

HUMAN PTGER2 PROTEIN, HFC TAG

货号: 11711

产品全名: 人 PTGER2 蛋白

规格: 10/50/100 µg

基因符号 EP2

目标蛋白: PTGER2

UNIPROT ID: P43116

描述: Recombinant Human PTGER2 Protein with C-terminal human Fc tag

背景: This gene encodes a receptor for prostaglandin E2, a metabolite of arachidonic acid which has different biologic activities in a wide range of tissues. Mutations in this gene are associated with aspirin-induced susceptibility to asthma. [provided by RefSeq, Oct 2009]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 28.7 kDa after removal of the signal peptide. The apparent molecular mass of PTGER2-hFc is approximately 35-55 kDa due to glycosylation.

分子特征: PTGER2(Met1-Glu23) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human PTGER2 Protein, hFc Tag on SDS-PAGE under reducing condition.