

HUMAN RSPO1(21-263) PROTEIN, HIS TAG

货号: 11567

产品全名: 人 RSPO1(21-263) 蛋白

规格: 10/50/100 µg

基因符号 R-spondin-1;hRspo1

目标蛋白: RSPO1

UNIPROT ID: Q2MKA7

描述: Recombinant human RSPO1(21-263) protein with C-terminal 6xHis tag

背景: This gene encodes a secreted activator protein with two cysteine-rich, furin-like domains and one thrombospondin type 1 domain. The encoded protein is a ligand for leucine-rich repeat-containing G-protein coupled receptors (LGR proteins) and positively regulates the Wnt signaling pathway. In mice, the protein induces the rapid onset of crypt cell proliferation and increases intestinal epithelial healing, providing a protective effect against chemotherapy-induced adverse effects. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 27.6 kDa after removal of the signal peptide. The apparent molecular mass of RSPO1(21-263)-His is approximately 35-55 kDa due to glycosylation.

分子特征: RSPO1(21-263)(Ser21-Ala263) 6×His tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human RSP01(21-263) Protein, His Tag on SDS-PAGE under reducing condition.