

## HUMAN SCGB2A2 PROTEIN, HFC TAG

货号: 11768

产品全名: 人 SCGB2A2 蛋白

规格: 10/50/100 µg

基因符号 MGB1;UGB2

目标蛋白: SCGB2A2

**UNIPROT ID:** Q13296

**描述:** Recombinant human SCGB2A2 protein with C-terminal human Fc tag

**背景:** Mammary gland specific. Over-expressed in breast cancer.

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 33.77 kDa after removal of the signal peptide.

**分子特征:** SCGB2A2 (Gly19-Phe93) hFc(Glu99-ALA330)

**纯化:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.