

HUMAN TMIGD2 PROTEIN, HFC TAG

货号: 11519

产品全名: 人 TMIGD2 蛋白

规格: 10/50/100 µg

基因符号 CD28 homolog;IGPR-1

目标蛋白: TMIGD2

UNIPROT ID: Q96BF3

描述: Recombinant human TMIGD2 protein with C-terminal human Fc tag

背景: Plays a role in cell-cell interaction, cell migration, and angiogenesis. Through interaction with HHLA2, costimulates T-cells in the context of TCR-mediated activation. Enhances T-cell proliferation and cytokine production via an AKT-dependent signaling cascade.[UniProtKB/Swiss-Prot Function]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 40.1 kDa after removal of the signal peptide. The apparent molecular mass of TMIGD2-hFc is approximately 55-70 kDa due to glycosylation.

分子特征: TMIGD2(Leu23-Gly150) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human TMIGD2 Protein, hFc Tag on SDS-PAGE under reducing condition.