

HUMAN TNFA PROTEIN, HIS TAG

货号: 11199

产品全名: 人 TNF α 蛋白

规格: 10/50/100 μ g

基因符号 DIF;NF- α ;NFA;NFSF2;achexin;achectin;NFA

目标蛋白: TNFA

UNIPROT ID: P01375

描述: Recombinant human TNF α protein with C-terminal 6xHis tag

背景: This gene encodes a multifunctional proinflammatory cytokine that belongs to the tumor necrosis factor (TNF) superfamily. This cytokine is mainly secreted by macrophages. It can bind to, and thus functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/TNFR. This cytokine is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation, apoptosis, lipid metabolism, and coagulation. This cytokine has been implicated in a variety of diseases, including autoimmune diseases, insulin resistance, and cancer. Knockout studies in mice also suggested the neuroprotective function of this cytokine.

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 18.1 kDa after removal of the signal peptide.

分子特征: TNF α (Val77-Leu233) 6xHis tag

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human TNF α Protein, His Tag on SDS-PAGE under reducing condition.

Human TNF α , His Tagged protein ELISA
0.1 μ g of Human TNF α , His Tagged protein per well




Figure 2. ELISA plate pre-coated by 1 μ g/ml (100 μ l/well) Human TNF α , His tagged protein (11199) can bind Anti-TNF α Neutralizing antibody 28064 in a linear range of 0.64-16 ng/ml.